



Contents lists available at ScienceDirect

Journal of Pharmaceutical Sciences

journal homepage: www.jpharmsci.org

Pharmaceutics, Drug Delivery and Pharmaceutical Technology

Storage stability and solution binding affinity of an Fc-fusion mimetic

Sama Prikalkhoran^{a,b}, David Guiliano^b, Hanieh Khalili^{a,c,*}^a School of Medicine and Biosciences, University of West London, W55RF, UK^b School of Life Sciences, University of Westminster, W1W 6UW, UK^c School of Pharmacy, University College London, WC1N 1AX, UK

ARTICLE INFO

Article history:

Received 7 September 2024

Revised 26 November 2024

Accepted 26 November 2024

Available online xxx

Keywords:

Fc-fusion mimetic

Aflibercept

MST

Binding affinity

Solution stability

Protein therapeutics

Bis-alkylation

ABSTRACT

This study evaluates the storage stability and solution binding affinity of a novel Fc-fusion mimetic, receptor-PEG-receptor (RpR), designed to address limitations of the current therapeutic aflibercept, a gold-standard therapy for age-macular degeneration (AMD). Using di(*bis*-sulfone) PEG linker as a structural scaffold, the mimetic aims to improve the storage stability and binding efficacy of the Fc fusion protein. Mass photometry and size-exclusion chromatography demonstrated that RpR, even in an unformulated buffer, exhibits superior storage stability exceeding 10 months compared to aflibercept. Furthermore, microscale thermophoresis was employed to determine RpR's binding affinity to VEGF in solution, providing a more physiologically relevant assessment than traditional binding assays. These findings highlight RpR's potential as a therapeutic candidate for the treatment of AMD disease, warranting further investigation.

Crown Copyright © 2024 Published by Elsevier Inc. on behalf of American Pharmacists Association. This is an open access article under the CC BY license (<http://creativecommons.org/licenses/by/4.0/>)

Introduction

The development of therapeutic proteins, particularly Fc-fusion proteins, has significantly advanced the treatment of many chronic diseases, including autoimmune disorders, cancers, and ocular diseases.^{1–3} Fc-fusion proteins combine the functional domain of a biologically active protein with the Fc region of an antibody.³ They offer therapeutic advantages through bivalency similar to IgGs, but they can be difficult to produce during early preclinical research and to scale for production.⁴ They are also prone to aggregation during downstream processing and have similar stability concerns as IgGs.^{5,6} In certain therapeutic contexts, such as organ-specific applications like ocular treatments, the Fc region may be unnecessary or even detrimental, particularly in managing inflammatory conditions.⁷ This highlights the need for alternative formats that retain the benefits of Fc-fusion proteins while addressing their limitations.

Aflibercept (Eylea) is a prominent example of an Fc-fusion protein that has revolutionised the treatment of neovascular age-related macular degeneration (AMD), diabetic macular edema, and other retinal vascular diseases. Aflibercept functions by binding to vascular endothelial growth factor (VEGF) and placental

growth factor (PlGF), preventing these factors from interacting with their receptors on the surface of endothelial cells and thereby inhibiting pathological angiogenesis and vascular permeability.⁸ Despite its success, aflibercept faces several challenges, including stability issues and the need for frequent intravitreal injections, which can lead to patient discomfort and increased risk of complications.^{9,10}

Structurally, aflibercept shares similarities with IgG antibodies due to the presence of an Fc domain. However, unlike IgG antibodies with their two heavy and two light chains, aflibercept is a homodimer glycoprotein consisting of two identical monomers. Each monomer is composed of vascular endothelial growth factor receptor 1 (VEGFR₁) and vascular endothelial growth factor receptor 2 (VEGFR₂) binding domains fused to the C2 and C3 regions of the Fc domain (VEGFR₁-VEGFR₂-Fc).⁸ These binding domains are linked by a disulfide bond, as shown in Fig. 1. The protein molecular weight of aflibercept is 96.9 kDa and contains 15 % glycosylation sites to give a total molecular weight of 115 kDa. It is a highly glycosylated protein, where 4 sites are located in the VEGFR domain and 1 site is located in the Fc domain.¹¹

To modify this structure, a di(*bis*-sulfone) PEG reagent can specifically target and react with the disulfide bond's cysteine thiols. This reaction replaces the disulfide bond with a stable 3-carbon methylene bridge. The PEG di(*bis*-sulfone) linker allows for site-specific

* Corresponding author.

E-mail address: hanieh.khalili@uwl.ac.uk (H. Khalili).




Fig. 1. Preparation of RpR from aflibercept using bis-alkylation bridging conjugation. The workflow shows the stages of RpR preparation beginning with the isolation of the VEGFR₁-VEGFR₂ dimer using proteolytic digestion and bis-alkylation bridging conjugation using PEG di(*bis*-sulfone) linker to produce RpR.

conjugation, providing structural stability while preserving the functional integrity of the protein.

Previously, we developed a novel Fc-fusion mimetic, termed Receptore-PEG-Receptor (RpR),¹² and antibody mimetic termed Fab-PEG-Fab (FpF),¹³ utilising a PEG di(*bis*-sulfone) linker as a scaffold to replace PEG with the Fc where two binding domains are linked together as if each binding domain is bound at the end of a linear molecule. FpF as IgG antibody mimetic and RpR as Fc-fusion mimetics are designed to have enhanced stability and binding properties compared to their parents' antibodies (Figs. 1, 2). Indeed, FpFs demonstrated superior protein stability compared to the parent IgG, with no aggregation or light and heavy chain dissociation observed in either liquid or lyophilised forms.¹⁴ However, the long term storage stability and binding properties of RpR in solution have not yet been investigated, which is the primary focus of this study.

In our earlier studies, we investigated the binding affinity of RpR using surface plasmon resonance (SPR) technology, where RpR's interaction with immobilised VEGF was analysed. These initial findings demonstrated promising binding characteristics and lower binding affinity for RpR compared with aflibercept.¹² However, SPR immobilises one binding partner, which does not fully replicate the dynamic interactions that occur in a physiological environment where both molecules are free in solution.

The primary objective of this current study is two-fold: firstly, to assess the stability of RpR under various storage conditions, including storage (4 °C) and physiological (37 °C) temperatures, and secondly, to evaluate its binding affinity to VEGF when both molecules

are in solution, closely mimicking the conditions within the human body. Stability is a critical factor for the efficacy and safety of therapeutic proteins, as instability can lead to aggregation, denaturation, and loss of function. These changes can reduce therapeutic effectiveness and increase the risk of adverse immune reactions. To address these objectives, we employed mass photometry to assess the stability of RpR over time at different temperature, tracking any potential aggregation or degradation. Mass photometry offered by Refeyn Ltd. is a relatively new technology that provides highly specific information about the mass and stoichiometry of biomolecules in solution. It is important to note that while mass photometry provides valuable insights into species distribution, it may over-represent smaller fragments due to their faster diffusion to the sensor. However, in this study, mass photometry served as a valuable tool for rapidly assessing the overall stability profiles of RpR and aflibercept at 37 °C. The primary goal is to compare their relative stability under different storage conditions, and mass photometry complemented by SDS-PAGE and size-exclusion chromatography (SEC) allowed us to track changes in the abundance of different species over time. Additionally, microscale thermophoresis (MST) was utilised to measure the binding affinity of RpR to VEGF in solution, providing a more physiologically relevant assessment compared to the previous SPR analysis.

The results indicate that RpR exhibits superior storage stability and solution binding affinity, positioning it as a promising candidate for further development in ocular inflammation therapies. These findings support the continued exploration of PEG di(*bis*-sulfone) conjugation as a valuable tool in the design of next-generation therapeutic proteins. In summary, By addressing the limitations of current treatments, such as aflibercept, we hope to pave the way for more effective and patient-friendly therapies for ocular and other chronic diseases.

Results and discussion

RpR preparation: Building upon previous work, aflibercept was subjected to proteolytic digestion using IdeS, cleaving below the hinge region while preserving the essential disulfide bridge for site-specific conjugation. This yielded a VEGFR₁-VEGFR₂ dimer (Fig. 3, lane 3), further purified using size exclusion chromatography (SEC) to prepare it for subsequent conjugation.

The site-specific conjugation process involved reducing the VEGFR₁-VEGFR₂ dimer to monomers using DTT as a reducing agent. Following removal of excess DTT via a PD-10 desalting column, the monomers were incubated with the PEG di(*bis*-sulfone) reagent (Fig. 3, lane 4). The resulting RpR was purified using SEC (Fig. 3, lane 5) and treated with sodium triacetoxyborohydride (STAB). This incubation aids in reducing the reactive ketone group to a more stable hydroxyl group,¹⁵ thereby minimising potential side reactions and improving the overall stability of the RpR conjugate.




Fig. 2. Structure of (A) IgG versus FpF, (B) Fc-fusion (aflibercept) versus RpR.




Fig. 3. Representative SDS-PAGE analysis for RpR preparation using di(*bis*-sulfone) PEG reagent. SDS-PAGE was stained with instant blue for protein staining. *lane 1*: Novex pre-stained protein marker, *lane 2*: Aflibercept (115 kDa), *lane 3*: VEGFR₁-VEGFR₂ dimer resulted from Ides digestion of aflibercept, *lane 4*: Reaction mixture between PEG di(*bis*-sulfone) reagent (1 eq) and reduced-VEGFR₁-VEGFR₂ dimer after 3 h incubation at ambient temperature, The slight difference in molecular weight between the two VEGFR₁-VEGFR₂ monomers might be attributed to variations in glycosylation.¹¹ *lane 5*: Purified RpR.

SDS-PAGE analysis confirmed that aflibercept migrated to an approximate molecular weight of 115 kDa (Fig. 3, lane 2). Incubation of aflibercept with IdeS enzyme (FabRICATOR, Genovis) specifically cleaved the Fc-fusion protein at the glycine-glycine bonds within the hinge region, yielding a VEGFR₁-VEGFR₂ dimer with a molecular weight of approximately 60 kDa. This was confirmed by SDS-PAGE analysis (Fig. 3, lane 3). Incubation of the VEGFR₁-VEGFR₂ dimer with DTT resulted in the reduction of the disulfide bond, yielding two lower molecular weight fragments of approximately 30 kDa each. These fragments are thought to be the desired VEGFR₁-VEGFR₂ monomer, each possessing a free thiol group available for subsequent bis-alkylation conjugation. The RpR conjugate was prepared by adding PEG di(*bis*-sulfone) reagent (1 equivalent) to the VEGFR₁-VEGFR₂ monomer and 3 h incubation at ambient temperature. Subsequent purification using SEC yielded the purified RpR, which migrated as a single band at approximately 70–80 kDa on SDS-PAGE (Fig. 3, lane 5).

RpR Stability Study: To assess the stability of RpR compared to aflibercept, both molecules were subjected to accelerated degradation studies by incubation at 37 °C for 1, 3, and 30 days. Mass photometry, a single-molecule imaging technique, was employed to monitor the structural integrity and potential aggregation of the proteins over time.^{16,17} This method allows for the direct visualisation and quantification of individual protein molecules based on their mass, providing valuable insights into their stability under physiological conditions.

Mass photometry analysis (Fig. 4, A) revealed a single peak for RpR at 78 kDa, confirming the purity and solution stability of the prepared conjugate. This result is consistent with the SDS-PAGE analysis (Fig. 3, lane 5), and previously published MALDI-ToF analysis which also showed a single band at approximately 70–80 kDa, further validating the successful purification and preparation of RpR. BAM (β -amylase) serves as a standard in mass photometry, offering multiple, distinct mass peaks that facilitate the calibration of the mass photometer's response.

Both RpR in an unformulated buffer (PBS only) and aflibercept (in its formulated buffer) exhibited stability at a concentration of 40 mg/mL when stored at 4 °C for a minimum of 60 days (Fig. 4, A). This finding aligns with previous reports on repackaged ziv-aflibercept, demonstrating maintained stability at 4 °C for up to 60 days.⁹ Upon

incubation at 37 °C, aflibercept exhibited early signs of instability, showing trace amounts of aggregation (high molecular weight peak) and chain dissociation even after 24 h (Fig. 4B). Conversely, RpR demonstrated remarkable stability, maintaining a single peak at 78 kDa for the entire 30-day incubation period at 37 °C. Long-term stability was studied using size-exclusion chromatography (SEC) and SDS-PAGE analysis. Both RpR (in PBS only) and aflibercept (formulation solution) were stored at 4 °C for 10 months and then analysed (Fig. 4, C). SEC and SDS-PAGE analysis revealed the presence of a high molecular weight peak (HMW) for aflibercept, indicative of aggregation. In contrast, RpR appeared as a single peak, demonstrating superior stability under these storage conditions.

Interestingly, RpR eluted earlier than aflibercept in SEC despite its lower molecular weight. It is important to note that SEC separation can be influenced by the specific column used, including the properties of the stationary phase. Variations in stationary pore size distribution and column material could lead to differences in elution profiles for the same molecules. AdvanceBio SEC column (2.7 μ m) was used to analyse the aggregation and/or degradation in RpR compared with aflibercept. This column has a unique hydrophilic coating to minimise secondary interactions between the sample and stationary phase.¹⁸ It is plausible that RpR exhibits minimal interaction with the column material, leading to faster elution, while aflibercept may experience some weak interactions that slightly delay its progress.

This enhanced stability of RpR is likely attributed to the replacement of the hinge region with the PEG di(*bis*-sulfone) linker. This flexible, 3-carbon bridge formed between cysteine residues, along with the formation of thio-ether bonds, imparts greater resilience to the molecule, preventing aggregation and chain dissociation. Similar stability improvement has been observed in PfF molecules compared to their parent IgG counterparts,¹⁴ suggesting that this modification strategy may be a broadly applicable approach for improving the stability of protein therapeutics.

RpR Binding Affinity in solution: Microscale thermophoresis (MST) was employed to evaluate the binding affinity of RpR and aflibercept to VEGF₁₆₅ in solution, offering a more physiologically relevant assessment compared to SPR techniques. MST's advantage lies in its ability to measure interactions in a free solution environment, eliminating potential artefacts caused by immobilisation and providing a more accurate representation of binding behaviour *in vivo*.^{19,20}



(A) Mass Photometry at 4°C**(B) Mass Photometry at 37°C**

Fig. 4. (A and B) Mass photometry analysis of aflibercept (115 kDa) and RpR (78 kDa) stability. BAM (β -amylase) is used as a standard which exists in solution as a mixture of monomers (58 kDa), dimers (109 kDa), and tetramers (225 kDa). (A) at 4 °C after 60 days, both aflibercept and RpR remain stable, showing no aggregation or chain dissociation. (B) at 37 °C after 30 days, RpR maintained its stability, while aflibercept displays high molecular weight (HMW) species and chain dissociation after 24 h. (C) Size-exclusion chromatography (SEC) and SDS-PAGE analysis to determine the stability of aflibercept (0.24 mg/mL in formulated buffer) and RpR (0.24 mg/mL in PBS only) after storage at 4 °C for 10 months. SEC analysis of aflibercept revealed a prominent high molecular weight (HMW) peak, indicative of aggregation, while RpR appeared as a single peak. SDS-PAGE analysis confirmed the presence of HMW species in the aflibercept sample, further demonstrating the superior stability of RpR under these storage conditions.

In this study, VEGF₁₆₅ was fluorescently labelled using an amine coupling reaction. While amine coupling is a widely used labelling method for MST, it's important to acknowledge its potential limitations. Heterogeneous labelling and steric hindrance due to the attached dye could influence the measured binding affinity.²¹

However, in this study, we chose to label VEGF to maintain consistency and minimise variability between the aflibercept and RpR binding measurements. VEGF. This approach also facilitated a direct head-to-head comparison between of the binding affinities of RpR and aflibercept to VEGF. Labelling either RpR or aflibercept would have

(C) Long-term stability study, 10 months at 4 °C




Fig. 4. Continued.

introduced potential bias due to the non-specific nature of amine coupling in the labelling process. This could potentially interfere with the binding sites, leading to inaccurate affinity measurements.

The degree of labeling (DOL) of VEGF₁₆₅ was optimised at 0.56 to ensure a strong and specific signal while minimising potential interference with the binding interaction. A range of concentrations of both RpR and aflibercept (0.9 μM to 0.00011 μM) was tested against a constant concentration of labelled VEGF₁₆₅ (0.625 μM) in MST capillaries to generate binding curves and accurately determine the binding affinity (KD) for each interaction. Fig. 5 demonstrates the binding curves and corresponding binding affinity for the interaction of both RpR and aflibercept with VEGF₁₆₅.

MST dose-response curves revealed that RpR had a lower dissociation constant ($K_D = 0.071 \text{ M}$) compared to aflibercept ($K_D = 0.230 \text{ M}$), indicating a higher binding affinity for VEGF₁₆₅. This observation is consistent with the SPR data,¹² which also demonstrated a slower dissociation rate (kd) for RpR compared to aflibercept, and higher binding affinity (K_D) indicating a longer residence time on the target. The enhanced affinity of RpR might be attributed to the replacement of the hinge region with a flexible PEG linker. This

modification may provide greater conformational freedom for the VEGFR₁-VEGFR₂ domains to interact optimally with VEGF₁₆₅, thus facilitating stronger and more stable binding. The flexibility of the PEG linker might also contribute to a faster association rate (ka) observed in the SPR data,^{12,13} allowing RpR to bind more rapidly to VEGF₁₆₅ compared to aflibercept.

The combination of MST and SPR data provides a comprehensive understanding of the binding kinetics and affinity of RpR and aflibercept towards VEGF₁₆₅. The results highlight the potential of modifying the linker region to enhance the binding properties of antibody-based drugs, offering a promising strategy for improving their therapeutic efficacy. The increased affinity of RpR, combined with its favourable stability profile, suggests that it may be a more potent and effective therapeutic option compared to aflibercept for the treatment of diseases mediated by VEGF₁₆₅. Previously published *in-vitro* angiogenesis assays of human umbilical vein endothelial cell (HUVEC) co-cultures indicated that RpR and aflibercept behave similarly in anti-angiogenic bioassays.¹² However, a more comprehensive set of *in vitro* and *in vivo* studies will be required to fully characterise RpR's anti-angiogenic and pharmacological properties.


MST measurement (0.9 μM to 0.00011 μM) at 37°C

Fig. 5. Microscale thermophoresis analysis comparing the binding affinities of aflibercept (red dots) and RpR (green dots) to VEGF₁₆₅. Serial dilutions of aflibercept and RpR (0.9 μM to 0.00011 μM) were tested against a constant concentration of fluorescently labelled VEGF₁₆₅ (0.625 μM). The resulting dose-response curves illustrate the binding behaviour and affinity (KD) of each molecule. Data points represent the mean of three independent experiments.

Conclusion

Building upon our previous successful generation and characterisation of RpR, this study further demonstrates its enhanced stability and binding affinity to VEGF. By replacing the Fc domain with a flexible PEG linker, we achieved several key advantages. First, we observed improved storage stability compared to the parent aflibercept, as evidenced by mass photometry analysis after incubation at 37 °C, and SEC and SDS-PAGE analysis of RpR stored in PBS buffer for over 10 months at 4 °C. This analysis revealed the remarkable long-term stability of RpR with no evidence of aggregation or PEG dissociation. Second, this modification reduces the potential for immunogenicity and avoids Fc-mediated effector functions. While these functions can be desirable in some therapeutic contexts (e.g., cancer immunotherapy), they are undesirable in the treatment of ocular diseases, where the eye is an immune-privileged site. Activating these effector functions in the eye could lead to inflammation and tissue damage. The substitution of the Fc domain with PEG in RpR aligns with the therapeutic goal of neutralising VEGF without triggering unwanted immune responses. This is particularly important considering the chronic nature of AMD and the need for repeated intravitreal injections.

In summary, the replacement of the Fc domain with a PEG scaffold in RpR offers a promising strategy for treating ocular diseases by mitigating immune-related side effects, potentially improving long-term safety and efficacy, and enhancing storage stability. This approach holds significant promise for addressing the limitations of current protein-based drugs and expanding the therapeutic landscape for a wide range of diseases.

Experimental section

Proteolytic digestion of aflibercept to prepare the dimeric VEGFR₁-VEGFR₂ fragment

Aflibercept (6.0 mg) was donated from Moorfield Eye Hospital after patient injections. It was digested using immobilised IdeS enzyme (FabRICATOR®, FragIT MidiSpin, Genovis; Cat no A0-FR6–100), following the optimised protocol reported previously in¹² and the manufacturer's instruction. Briefly, the column was equilibrated with cleavage buffer (50 mM sodium phosphate, 150 mM NaCl, pH 6.6) and aflibercept (6 mg in 1.0 mL cleavage buffer) was loaded onto the column. The digestion was incubated for 30 min at room temperature with end-over-end mixing. The resulting digestion mixture was then purified using a CaptureSelect MidiSpin column (Genovis) packed with a multi-species Fc affinity matrix, according to the manufacturer's instructions. The purified VEGFR₁-VEGFR₂ dimer was eluted, analysed by SDS-PAGE, and quantified using a micro BCA assay, yielding 2 mg of purified protein.

Preparation of RpR general procedure

Dimeric VEGFR₁-VEGFR₂ fragment (0.8 mg) was reduced with dithiothreitol (DTT, 6.0 mM) in PBS (pH 7.3) for 30 min at room temperature to yield monomeric VEGFR₁-VEGFR₂. The reaction mixture was buffer exchanged into sodium phosphate buffer (20 mM), EDTA (10 mM), pH 7.6 using a PD-10 column to remove excess DTT. The monomeric VEGFR₁-VEGFR₂ (0.24 mg/mL, 3.3 mL) was then conjugated with 0.9 equivalents of a 10 kDa PEG di(*bis*-sulfone) reagent (previously described^{12,13}) for 12 h at 4 °C. The resulting RpR was purified by size exclusion chromatography (Superose 12 HR 10/30 column, PBS mobile phase) and fractions containing the desired product were pooled, concentrated, and quantified by micro BCA assay, yielding 0.2 mg of purified RpR.

Determination of the mass size of RpR by mass photometry

The mass photometry (Refeyn Ltd) was used to monitor molecular weight, aggregation, and degradation of the protein over time. Samples (aflibercept and RpR) were prepared at a concentration of 0.3 mg/mL using PBS buffer, which was passed through a 0.2 μm filter. Briefly, samples and standards (β-Amylase and BSA, both 100 nM) were diluted in PBS and analyzed using a mass photometer. The instrument was prepared by applying immersion oil to the objective lens and positioning magnets on the stage. Samples were introduced using the droplet dilution method (2 μL sample in 18 μL PBS) and analysed for 60 s, during which molecular weight and particle counts were recorded. Stability samples were analysed at specified time intervals following incubation.

Determination of stability using Size-Exclusion Chromatography and SDS-PAGE analysis Size-exclusion chromatography (SEC) was conducted to evaluate the stability of RpR and aflibercept, specifically to detect any degradation or chain dissociation occurring over 10 months of storage at 4 °C. The samples included RpR at a concentration of 0.24 mg/mL and aflibercept, diluted to a final concentration of 0.24 mg/mL in PBS. Samples were injected into an Agilent 1260 Infinity HPLC system, using PBS (pH 7.3) as a mobile phase and a flow rate of 0.35 mL/min. Each run was completed over 35 min, and separation was achieved using an AdvanceBio SEC column (300 Å, 2.7 μm, 7.8 × 300 mm). Detection was monitored at 280 nm to accurately capture any potential structural changes or dissociations in the protein samples.

Non-reducing SDS-PAGE was employed to analyse protein samples without disrupting disulfide bonds. Samples (20 μL) were mixed with a non-reducing sample buffer (4 μL, Pierce™ LDS Sample Buffer, Non-Reducing (4X)), and loaded (10 μL) loaded into precast SDS-PAGE gels (Novex™ Tris-Glycine Mini Protein Gels, 4–20 %, 1.0 mm). Electrophoresis was conducted at a constant voltage of 160 V for 1 hour using NuPAGE™ MOPS SDS Running Buffer (20X). The gels were stained with Coomassie blue for protein staining and a silver staining kit (Pierce™ Silver Stain Kit) to visualise any trace quantity of protein.

Determination of the solution binding affinity of RpR by microscale thermophoresis

Microscale thermophoresis (Monolith, NanoTemper) was used to determine the binding affinity of aflibercept, and RpR to human recombinant VEGF₁₆₅ (38 kDa). VEGF₁₆₅ (10 μM) was labelled with a fluorescent dye according to the manufacturer's instructions (Protein Labeling Kit RED-NHS 2nd Generation (Amine Reactive), Nanotemper) and purified using the provided purification column (Monolith premium capillaries, included in the kit). Briefly, VEGF₁₆₅ (10 μM in labelling buffer) was incubated with the RED-NHS dye (300 μM in labelling buffer containing 50 % DMSO) for 30 min at room temperature in the dark. The labelled protein was purified using the provided gravity flow column (Monolith premium capillaries, NanoTemper Technologies) with phosphate-buffered saline (PBS) as the elution buffer. The degree of labelling (DOL), defined as the ratio of dye molecules to protein molecules, was determined spectrophotometrically by measuring the absorbance at 205 nm (protein) and 650 nm (dye) and was optimised to fall within the range of 0.5–1.0 for optimal MST signal. A systematic approach was used to optimise the DOL by varying the dye-to-protein ratio in a series of small-scale labelling reactions and assessing the DOL after each purification.

For the MST assay, serially diluted RpR and aflibercept samples were prepared in PBS (0.9 μM to 0.00011 μM) and mixed with labelled VEGF₁₆₅ (0.15 μM) in glass capillaries. MST measurements were performed using a Monolith Pico system with MST power of 40 % and excitation power of 80 %. The dissociation binding constant

(KD) for each interaction was determined using the MST analysis software.

Declaration of competing interest

The authors declare that they have no known competing financial interests or personal relationships that could have appeared to influence the work reported in this paper.

Acknowledgements

S.P. is thankful for the Vice Chancellor PhD scholarship from the University of West London. S.P. is also grateful for having access to mass photometry at UCL School of Pharmacy and UCL Institute of Neurology. *Special thanks to Dr. Maria Kalyukina, Department of Clinical and Experimental Epilepsy, to assist with mass photometry analysis.

References

1. Ultee ME. Therapeutic Fc-Fusion proteins: edited by Steven M Chamow, Thomas Ryll, Henry B Lowman and Deborah Farson. *MAbs*. 2014;6(4):810–811. <https://doi.org/10.4161/mabs.28832>. Epub. 2014 Apr 8. PMID: PMC4171016.
2. Duivelshof BL, Murisier A, Camperi J, et al. Therapeutic Fc-fusion proteins: current analytical strategies. *J Sep Sci*. 2021;44(1):35–62. <https://doi.org/10.1002/jssc.202000765>. Epub. 2020 Sep 29. PMID: 32914936.
3. Strohl WR. Fusion proteins for half-life extension of biologics as a strategy to make biobetters. *BioDrugs*. 2015;29(4):215–239. <https://doi.org/10.1007/s40259-015-0133-6>. PMID: 26177629; PMCID: PMC4562006.
4. Hoogenboom HR. Selecting and screening recombinant antibody libraries. *Nat Biotechnol*. 2005;23(9):1105–1116. <https://doi.org/10.1038/nbt1126>. PMID: 16151404.
5. Wang W. Advanced protein formulations. *Protein Sci*. 2015;24(7):1031–1039. <https://doi.org/10.1002/pro.2684>. Epub. 2015 May 1. PMID: 25858529; PMCID: PMC4500304.
6. Cordes AA, Platt CW, Carpenter JF, Randolph TW. Physical instability of a therapeutic Fc fusion protein: domain contributions to conformational and colloidal stability. *J Pharm Sci*. 2009;48(49):11724–11736. <https://doi.org/10.1021/bi900853v>. Biochemistry.
7. van der Neut Kofschoten M, Schuurman J, Losen M, et al. Anti-inflammatory activity of human IgG4 antibodies by dynamic Fab arm exchange. *Science*. 2007;317(5844):1554–1557. <https://doi.org/10.1126/science.1144603>. PMID: 17872445.
8. EMA (European Medicines Agency). Assessment Report For Aflibercept (EMA/646256/2012), European Medicines Agency Committee for Medicinal Products for Human Use, Editor., 2012, p. 83. Available from: <https://www.ema.europa.eu/en/homepage>.
9. Juhong J, Pongsachareonnon PF, Somkijrungraj T, et al. The sterility, stability and efficacy of repackaged ziv-aflibercept for intravitreal administration. *Sci Rep*. 2022;12(1):2971. <https://doi.org/10.1038/s41598-022-06831-2>. PMID: 35194061; PMCID: PMC8863784.
10. Moreno MR, Tabitha TS, Nirmal J, et al. Study of stability and biophysical characterization of ranibizumab and aflibercept. *Eur J Pharm Biopharm*. 2016;108:156–167. <https://doi.org/10.1016/j.ejpb.2016.09.003>. Epub 2016 Sep 8. PMID: 27615995.
11. Shen Z, Wang Y, Xu H, et al. Analytical comparability assessment on glycosylation of ziv-aflibercept and the biosimilar candidate. *Int J Biol Macromol*. 2021;180:494–509. <https://doi.org/10.1016/j.ijbiomac.2021.03.020>. Epub. 2021 Mar 6. PMID: 33684428.
12. Khalili H, Khaw PT, Brocchini S. Fc-fusion mimetics. *Biomater Sci*. 2016;4(6):943–947. <https://doi.org/10.1039/c6bm00077k>. Epub 2016 Apr 29. PMID: 27127811.
13. Khalili H, Godwin A, Choi JW, Lever R, Khaw PT, Brocchini S. Fab-PEG-Fab as a potential antibody mimetic. *Bioconjug Chem*. 2013;24(11):1870–1882. <https://doi.org/10.1021/bc400246z>. Epub. 2013 Nov 4. PMID: 24073593.
14. Khalili H. Using different proteolytic enzymes to digest antibody and its impact on stability of antibody mimetics. *J Immunol Methods*. 2021;489:112933. <https://doi.org/10.1016/j.jim.2020.112933>. Epub. 2020 Nov 21. PMID: 33232747.
15. AF Abdel-Magid, Carson KG, Harris BD, Maryanoff CA, Shah RD. Reductive amination of aldehydes and ketones with sodium triacetoxyborohydride. Studies on direct and indirect reductive amination procedures (1). *J Org Chem*. 1996;61(11):3849–3862. <https://doi.org/10.1021/jo960057x>. PMID: 11667239.
16. Sonn-Segev A, Belacic K, Bodrug T, et al. Quantifying the heterogeneity of macromolecular machines by mass photometry. *Nat Commun*. 2020;11(1):1772. <https://doi.org/10.1038/s41467-020-15642-w>. PMID: 32286308; PMCID: PMC7156492.
17. Wu D, Xu P, Kelly M, Ryan ET, Kováč P, Piszczek G. Mass photometry: a powerful tool for carbohydrates-proteins conjugation monitoring and glycoconjugates molecular mass determination. *Glycoconj J*. 2023;40(4):401–412. <https://doi.org/10.1007/s10719-023-10126-7>. Epub. 2023 Jul 1. PMID: 37392327; PMCID: PMC10374364.
18. <https://www.agilent.com/en/product/biopharma-hplc-analysis/aggregate-fragment-analysis/advancebio-sec>.
19. Jerabek-Willemsen M, Wienken CJ, Braun D, Baaske P, Dühr S. Molecular interaction studies using microscale thermophoresis. *Assay Drug Dev Technol*. 2011;9(4):342–353. <https://doi.org/10.1089/adt.2011.0380>. PMID: 21812660; PMCID: PMC3148787.
20. Wienken CJ, Baaske P, Rothbauer U, Braun D, Dühr S. Protein-binding assays in biological liquids using microscale thermophoresis. *Nat Commun*. 2010;1:100. <https://doi.org/10.1038/ncomms1093>. PMID: 20981028.
21. Jiang H, Cole PA. N-terminal protein labeling with N-hydroxysuccinimide esters and microscale thermophoresis measurements of protein-protein interactions using labeled protein. *Curr Protoc*. 2021;1(1):e14. <https://doi.org/10.1002/cpz1.14>. Erratum in: *Curr Protoc*. 2022 Aug;2(8):e552. 10.1002/cpz1.1552. PMID: 33484499; PMCID: PMC7839251.